

Electrophoresis and Sequencing Basics

Review: Mitochondrial lab




DNA Extraction from buccal cheek cells

Review: Mitochondrial DNA Structure



Human Mitochondrial DNA 16569 bp




Review: Mitochondrial lab




DNA Amplification by PCR

Review: Mitochondrial lab



Electrophoresis



Electrophoresis – movement by electricity


- Allows us to visualize DNA (PCR product)
- We can estimate the size (distance a "band" travels)

We can estimate the amount (brightness of the band


Electrophoresis – Components


- Gel tank: Holds electrolytic buffer and has electrodes where electricity is applied
- Buffer: A salt solution that conducts electricity
- Gel: Made of agar; holds the DNA and separates DNA according to size
- Loading Dye: Makes it easier to load the PCR sample and keeps the DNA in the gel.
- UV Dye: Binds to DNA. Can be in the gel or the loading dye. (We use Syber Green.
 Ethidium Bromide can also be used)
- We can estimate the amount (brightness of the band


Electrophoresis – Setup





Electrophoresis – Setup























Electrophoresis – Cut DNA example




Mitochondrial DNA – replicated by PCR





DNA Sequencing



DNA Sequencing

- Allows us to determine the nucleotides and their order on a segment of DNA
- There are many types of DNA sequencing. For PCR products we use "Sanger" (cycle) sequencing (Fredrick Sanger)
- In this method we usually can sequence up to 1000bp of DNA
- Many components in common with PCR so be careful not to confuse them!

DNA Sequencing

Cycle Sequencing

